

Hit and Miss Transform in Medical Image Processing.



Faculty-Dr. Tapan Kumar Jain

Report By-Divit Srivastava (BT23CSE224) Erukulla

Raghunandan (BT23CSE182) Arnab Halder

(BT23CSE034) Prathmesh Ghormare (BT23CSE039)

Atharva Khaire (BT23CSE042)

Github Link-(<https://github.com/bt23cse224/Hit-or-Miss-Transform-for-Tumor-Detection>)

Introduction

A basic morphological operation for shape detection, pattern matching, and structural analysis in binary images is the Hit-and-Miss Transform (HMT). It is very effective at identifying precise geometrical structures because it is the only basic morphological operation that can simultaneously detect particular foreground and background pixel configurations. Conceptually, HMT identifies whether a specific shape exists at a given location by analysing an image using two complementary structuring elements: one that must fit the foreground ("hit") and another that must fit the background ("miss").

Applications of HMT Transform in Brain Tumor Detection.


In brain tumor detection, the Hit-and-Miss Transform is used after generating a binary mask of the brain or tumor candidate regions. Tumors often exhibit irregular boundaries, disorganized tissue structures, and specific morphological discontinuities that differ significantly from normal brain anatomy. By designing structuring elements that capture these characteristics such as small protrusions, boundary irregularities, or texture-based patterns, the HMT can highlight subtle structural abnormalities. It helps refine tumor segmentation by detecting tumor edges, eliminating false positives, and identifying shape-specific tumor markers. In MRI-based tumor analysis, the HMT acts as a precise tool for validating morphological consistency and improving the accuracy of region-of-interest.

Objective:

The primary objective of this project is:

- To implement the Hit-or-Miss Transform for detecting tumor-like structures in MRI brain images.
- To design appropriate structuring elements (kernels) that approximate the shape of a tumor.
- To visualize and analyze the intermediate steps including binary conversion, erosions, and final detection.
- To produce a system that outputs whether a tumor-like pattern exists in the given MRI image.

Workflow



Detailed Process Explanation:

1. Image Loading and Grayscale Conversion

```
image_gray = cv2.imread(image_path, cv2.IMREAD_GRAYSCALE)
```

2. Binary Conversion

The grayscale image is thresholded to a binary form:

```
_, binary_image = cv2.threshold(image_gray, 150, 255,  
cv2.THRESH_BINARY)
```

- Pixel intensity $\geq 150 \rightarrow$ white (foreground)
- Pixel intensity $< 150 \rightarrow$ black (background)

This isolates bright tumor regions.

3. Structuring Element Design

3.1 Hit Kernel (Foreground)

A 20×20 elliptical kernel:

```
hit_kernel = cv2.getStructuringElement(cv2.MORPH_ELLIPSE, (20, 20))
```

This models tumor-like shapes which are typically round/elliptical.

3.2 Miss Kernel (Background)

A 40×40 elliptical kernel with its center carved out:

```
miss_kernel = cv2.getStructuringElement(cv2.MORPH_ELLIPSE, (40, 40))  
miss_kernel[10:30, 10:30] = 0
```

This ensures:

- The center matches the tumor region
- The surrounding area must be background

This combination enforces precise tumor detection.

4. Hit Pass – Erosion on Binary Image

```
erosion_hit = cv2.erode(binary_image, hit_kernel, iterations=1)
```

This identifies regions where the tumor-like shape fits perfectly inside the bright foreground.

5. Miss Pass – Erosion on Complement Image

The binary image is inverted:

```
complement_image = cv2.bitwise_not(binary_image)
```

Erosion with miss kernel:

```
erosion_miss = cv2.erode(complement_image, miss_kernel, iterations=1)
```

This identifies locations where the background pattern fits.

6. Combining Results – Hit-or-Miss Operation

```
hit_or_miss_result = cv2.bitwise_and(erosion_hit, erosion_miss)
```

Only those pixels where both conditions are satisfied are retained.

7. Visualization

The program displays five stages:

1. Original MRI
2. Binary Image
3. Hit Erosion
4. Miss Erosion

5. Final Hit-or-Miss Result



8. Final Detection

```
if np.any(hit_or_miss_result):  
    print("DETECTION: A tumor-like shape was found.")  
else:  
    print("NO DETECTION: No tumor-like shape found.")
```

Experimental Output

The output images showcase:

- The binary segmentation
- Regions matching tumor-like shapes (Hit result)
- Regions matching the background pattern (Miss result)



Advantages of This Approach:

- **Shape-based detection:** Useful when analyzing morphological abnormalities.
- **Low computational complexity:** Erosion and logical operations are fast.

-
- **Easy to interpret:** Structuring elements give explicit shape detection capability.

Conclusion:

The Hit-or-Miss Transform is a powerful morphological tool for detecting specific patterns in medical images. In this project, it was successfully applied to detect tumor-like shapes in brain MRI scans by designing appropriate structuring elements and analyzing both foreground and background patterns.

Future Scope

- Adaptive structuring elements based on image features.
- Integration with CNN-based tumor segmentation models.
- Automating threshold selection using Otsu's method.
- Applying morphological thinning before HMT to enhance shape consistency.